

Vagiwell®

Dilator 5 LS

Vaginal dilator for personal use

Indications

- neovagina (sex reassignment surgery), post-surgery
- lichen sclerosus (to keep the vaginal entrance functional)

Application

The Vagiwell® Dilator 5LS is cleaned with hot water and mild soap before first use and dried completely.

The entrance of the vagina and the dilator are moistened with a water-soluble lubricant and the dilator is inserted in a relaxed position.

For removal, the dilator is gently pulled out of the vagina.

Lichen sclerosus patients usually use the dilator once daily to keep the sclerotic (firm) entrance of the vagina and the vulva dilatable.

If desired, Vagiwell® Dilator sizes 1 – 5 can initially be used progressively. If Vagiwell® 5 LS is not sufficient, Vagiwell® No.6 is available as an alternative.

Cleaning

The Vagiwell® Dilator 5 LS is intended for personal use by one user only. After use, the dilator is cleaned with warm water and mild soap.

Content

- one Vagiwell® Dilator 5 LS
- water-soluble lubricant, 50 ml
- storage bag
- instruction manual

Material

- Vagiwell® Dilator: silicone
- Storage bag: cotton/nylon®

Vagiwell® Dilators 5 LS

Art. Nr.: Vagiwell 5 LS

EAN: 4013273002165

PZN: 15994727

Size	Colour	Diameter A / B	Length C
5 LS	white	33/36 mm	181 mm

