

Vagiwell®

Dilator 6

Vaginal dilator for personal use

Indications

- neovagina (sex reassignment surgery), post-surgery
- lichen sclerosus (in case that the Vagiwell® Dilator 5 LS is not sufficient)

Application

Prior to first use, the dilator should be cleaned with warm water and mild soap and then thoroughly dried. The vaginal opening and the dilator should be moistened with a water-soluble lubricant. The dilator should then be carefully inserted in a relaxed position.

To remove it, the dilator is carefully pulled out of the vagina.

In the event of pain, the dilator should be removed immediately and further use clarified with the treating physician.

If desired, the Vagiwell® Dilator sizes 1 – 5 and 5 LS can first be used in a progressively increasing manner before Vagiwell® 6 is used for the first time.

Cleaning

The Vagiwell® Dilator 6 is intended for personal use by one user only. After use, the dilator is cleaned with warm water and mild soap.

Content

- one Vagiwell® Dilator 6
- water-soluble lubricant, 50 ml
- storage bag
- instruction manual

Material

- Vagiwell® Dilator: silicone
- Storage bag: cotton/nylon®

Vagiwell® Dilator 6

Art. Nr.: Vagiwell 6

EAN: 4013273001632

PZN: 11008157

Size	Colour	Diameter A / B	Length C
No. 6	white	39/41 mm	188 mm

